



Anti-C1ORF177 (aa 252-301 9SMQKKKPRELMNFKSFVEELNSHHNKKHGVFSKLP TPTERIYWANL) polyclonal antibody (CABT- BL651)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within internal sequence amino acids 252-301 9SMQKKKPRELMNFKSFVEELNSHHNKKHGVFSKLP (TPTERIYWANL) of Human C1orf177
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [163747](#)

Protein Refseq [NP_001104003](#)

UniProt ID [Q3ZCV2](#)
