



# Anti-C1ORF38 (full length) polyclonal antibody (CABT-BL658)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Full length protein, corresponding to amino acids 1-132 of Human C1orf38 (AAH81568.1).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	1X PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">9473</a>
<b>Protein Refseq</b>	<a href="#">NP_001099026</a>
<b>UniProt ID</b>	<a href="#">Q5TEJ8</a>