



Anti-C5ORF32 monoclonal antibody, clone 5F22 (CABT-BL748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ORF1-FL49 (AAH13643, 1 a.a. ~ 98 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5F22
Conjugate	Unconjugated
Applications	WB,ELISA
Buffer	In 1x PBS, pH 7.2
Preservative	None
Storage	Store at -2°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	84418
Protein Refseq	NP_115788

UniProt ID

[Q9H1C7](#)
