



# Anti-25-OH Vitamin D3 monoclonal antibody, clone SN4 (CABT-BL8456)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Reacts to 25-OH Vitamin D3. It does not cross-react with Vitamin D3, Vitamin D2, and 25-OH Vitamin D2.
<b>Immunogen</b>	25-OH Vitamin D3 BSA conjugate
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Clone</b>	SN4
<b>Purification</b>	Protein A affinity purified from an animal origin-free culture supernatant.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IA
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Stable for 1 Year at -20.0°C from date of receipt