



Anti-CCDC3 (C-terminal) polyclonal antibody (CABT-BL862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide selected from the C-terminal region of Human Coiled Coil Domain Containing Protein 3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Endoplasmic reticulum
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Coiled-coil domain-containing protein 3 is a protein that in humans is encoded by the CCDC3 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [83643](#)

Protein Refseq [NP_113643](#)

UniProt ID [Q9BQI4](#)
