



Anti-Chikungunya E1 monoclonal antibody (CABT-BL8724)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No reactivity with Dengue and Zika viruses
Target	Chikungunya
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	CHIKV
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	0.85% Sodium Chloride
Preservative	None
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.