



Mouse Anti-Cryptosporidium monoclonal antibody, clone 8E1 (CABT-BL8752)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Cryptosporidium
Source/Host	Mouse
Species Reactivity	Cryptosporidium
Clone	8E1
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Lyophilized
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year