



Anti-6xHIS Tag monoclonal antibody (CABT-BL8767)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	6xHIS tag
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IF, WB, IHC Recommended dilution: ELISA: Assay Dependent IF: 10 - 2 ug/ml IHC: 5.0 - 1.0 ug/ml WB: 0.2 - 0.05 ug/ml
Format	Lyophilized
Size	100 µg
Buffer	Phospate buffered saline, pH 7.2
Preservative	None
Storage	4°C for 1 year