



Mouse Anti-SARS-CoV-2 Spike Neutralizing Monoclonal antibody, clone NN54 (CABT-CS064)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	React with SARS-CoV-2 Spike RBD Protein. No cross-reactivity in ELISA with SARS-CoV, MERS-CoV, HCoV-HKU1, HCoV-NL63, HCoV-229E, HCoV-OC43 Spike S1 Protein.
Target	SARS-CoV-2 Spike S1
Immunogen	Recombinant SARS-CoV-2 Spike S1-mFc Protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	SARS-CoV-2
Clone	NN54
Purification	Protein A
Conjugate	unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	None

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

BACKGROUND

Introduction

The spike (S) glycoprotein of coronaviruses contains protrusions that will only bind to certain receptors on the host cell. Known receptors bind S1 are ACE2, angiotensin-converting enzyme 2; DPP4, dipeptidyl peptidase-4; APN, aminopeptidase N; CEACAM, carcinoembryonic antigen-related cell adhesion molecule 1; Sia, sialic acid; O-ac Sia, O-acetylated sialic acid. The spike is essential for both host specificity and viral infectivity. The term 'peplomer' is typically used to refer to a grouping of heterologous proteins on the virus surface that function together. The spike (S) glycoprotein of coronaviruses is known to be essential in the binding of the virus to the host cell at the advent of the infection process.

Keywords

SARS-CoV-2; coronavirus; SARS-CoV-2 spike 1; SARS-CoV-2 spike protein; SARS-CoV-2 S1; SARS-CoV-2 RBD; SARS-CoV-2 spike RBD
