



Rabbit Anti-Human SBSN Polyclonal Antibody (CABT-L2220)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Polyclonal Antibody to Suprabasin (Knockout Validated)
Specificity	The antibody is a rabbit polyclonal antibody raised against SBSN. It has been selected for its ability to recognize SBSN in immunohistochemical staining and western blotting.
Target	SBSN
Immunogen	Recombinant fragment corresponding to human SBSN (Gly434~Pro590)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Pig
Purification	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	Lot specific
Size	200 µg
Buffer	Supplied as solution form in 0.01M PBS with 50% glycerol, pH7.4.
Preservative	0.05% Proclin-300

Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
Ship	4°C with ice bags

BACKGROUND

Introduction	SBSN (suprabasin) is a protein-coding gene.
Keywords	SBSN

GENE INFORMATION

Gene Name	SBSN suprabasin [Homo sapiens (human)]
Official Symbol	SBSN
Synonyms	SBSN; suprabasin; UNQ698; HLAR698;
Entrez Gene ID	374897
Protein Refseq	NP_001159506
UniProt ID	Q6UWP8
Chromosome Location	19q13.13
Function	molecular_function;