



Mouse Anti-Chromium monoclonal antibody, clone Cr (CABT-L2516)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Chromium
Immunogen	Chromium with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	Cr
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Chromium is a chemical element with symbol Cr and atomic number 24. It is the first element in group 6. It is a steely-grey, lustrous, hard and brittle metal which takes a high polish, resists tarnishing, and has a high melting point.

Keywords

Chromium;Cr
