



Mouse Anti-Citrinin monoclonal antibody, clone CIT (CABT-L2517)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Citrinin
Immunogen	Citrinin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CIT
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Citrinin is a mycotoxin which is often found in food. It is a secondary metabolite produced by fungi that contaminate long-stored food and it causes different toxic effects, like nephrotoxic, hepatotoxic and cytotoxic effects. Citrinin is mainly found in stored grains, but sometimes also in fruits and other plant products.

Keywords

Citrinin;
