



Mouse Anti-Medroxyprogesterone monoclonal antibody, clone MP (CABT-L2521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Medroxyprogesterone
Immunogen	Medroxyprogesterone with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MP
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Medroxyprogesterone (MP), is a progestin which is not used medically. A derivative, medroxyprogesterone acetate (MPA), is used as a medication in humans, and is far more widely known in comparison. Medroxyprogesterone is sometimes used as a synonym for medroxyprogesterone acetate, and what is almost always being referred to when the term is used is MPA and not medroxyprogesterone.

Keywords

Medroxyprogesterone;17 α -hydroxy-6 α -methylprogesterone;MP