



Mouse Anti-Spectinomycin monoclonal antibody, clone SPM (CABT-L2528)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Spectinomycin
Immunogen	Spectinomycin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	SPM
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Spectinomycin was discovered in 1961. It is on the World Health Organization's List of Essential Medicines, the most effective and safe medicines needed in a health system. It is made from the bacterium *Streptomyces spectabilis*.

Keywords

Trobicin; Spectinomycin; (2R;4aR;5aR;6S;7S;8R;9S;9aR;10aS)-4a;7;9-trihydroxy-2-methyl-6;8-bis(methylamino)decahydro-4H-pyrano[2;3-b][1;4]benzodioxin-4-one