



Mouse Anti-Paromomycin monoclonal antibody, clone PAM (CABT-L2529)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Cross reactivity: Neomycin 20ppb
Immunogen	Paromomycin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PAM
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

Paromomycin is an antibiotic used to treat a number of infections including amebiasis, giardiasis, leishmaniasis, and tapeworm infection. It is a first line treatment for amebiasis or giardiasis during pregnancy. Otherwise it is generally a second line treatment option. It is used by mouth, applied to the skin, or by injection into a muscle

Keywords

Catenulin;Aminosidine;monomycin;aminosidin
