



# Mouse Anti-Maduramicin monoclonal antibody, clone MAD (CABT-L2535)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Maduramicin
Immunogen	Maduramicin with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MAD
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

Maduramicin is a monoglycoside polyether derived from the fungus *Actinomadura rubra*.

---

**Keywords**

Maduramicin; Maduramycin; Maduramicin Ammonium; Maduramycin Ammonium; Cygro;

---