



Mouse Anti-4-Chlorophenoxyacetic acid monoclonal antibody, clone 4-D (CABT-L2543)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	4-Chlorophenoxyacetic acid
Immunogen	4-Chlorophenoxyacetic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	4-D
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	-20°C. Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

4-Chlorophenoxyacetic acid or parachlorophenoxyacetate (pCPA) is a synthetic pesticide similar to chemicals in a group of plant hormones called auxins.

Keywords

4-Chlorophenoxyacetic acid;2-(4-Chlorophenoxy)acetic acid;Parachlorophenoxyacetate;4-Chlorophenoxyacetate
