



Mouse Anti-Diaveridine monoclonal antibody, clone DVD (CABT-L3055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Diaveridine with carrier protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DVD
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

In some countries, this medicine may only be approved for veterinary use.

Keywords

Diaveridinum (Latin);Diaveridin (German);Diavéridine (French);Diaveridina (Spanish)