



Mouse Anti-Ormetoprim monoclonal antibody, clone OMP (CABT-L3056)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Immunogen | Ormetoprim with carrier protein. |
| Isotype | IgG |
| Source/Host | Mouse |
| Species Reactivity | N/A |
| Clone | OMP |
| Purification | Purified from mouse ascites. |
| Conjugate | Unconjugated |
| Applications | ELISA, LFIA |
| Format | Liquid |
| Concentration | Lot specific |
| Size | 1 mg |
| Buffer | PBS |
| Preservative | None |
| Storage | Long time storage is recommended at -20°C. |
| Ship | Wet ice |

BACKGROUND

Introduction

Ormetoprim is a close analog of diaveridine and a folic acid antagonist used in potentiated sulfonamide mixtures used as coccidiostats.

Keywords

Ormetoprim
