



Mouse Anti-Azithromycin monoclonal antibody, clone ATM (CABT-L3058)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Azithromycin with carrier protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ATM
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Azithromycin is an antibiotic useful for the treatment of a number of bacterial infections. This includes middle ear infections, strep throat, pneumonia, traveler's diarrhea, and certain other intestinal infections. It may also be used for a number of sexually transmitted infections including chlamydia and gonorrhea infections. Along with other medications, it may also be used for malaria

Keywords

9-deoxy-9 α -aza-9 α -methyl-9 α -homoerythromycin A;Azithromycin