



Chimeric Human Anti-PEG monoclonal antibody, clone Dbhq5 (CABT-L3141)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Fusion protein of mouse BHQ5 VH and VL with human IgM
Specificity	This is a chimeric antibody composed of the BHQ5 anti-PEG variable region genes fused to human IgM constant domains. It is expressed with the human J-chain to form pentamers.
Target	This antibody binds to the backbone of PEG and can be detected with anti-human secondary antibodies.
Isotype	IgM
Source/Host	Human
Species Reactivity	N/A
Clone	dBHQ5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	50 µg
Buffer	Purified antibody supplied in 50% glycerol and phosphate buffered saline (0.14 M NaCl, 2.7 mM KCl, 1.5 mM KH ₂ PO ₄ , 8.1 mM Na ₂ HPO ₄)

Preservative	0.04% Sodium Azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	<p>PEG (polyethylene glycol) is a water-soluble, nontoxic, biocompatible polymer that has been approved by the Food and Drug Administration (FDA) for human intravenous, oral and dermal applications. Attachment of PEG chains to proteins can reduce their immunogenicity, minimize proteolytic cleavage and increase their serum half-life. PEG has also been attached to small molecules and liposomes for more selective delivery. PEG-modification of superparamagnetic iron oxide and quantum dots can improve their biocompatibility and reduce non-specific uptake. PEG antibodies can be a vital tool for propelling therapeutics to market by serving as a positive control anti-drug antibody, measuring clearance of a drug, or simply as a QA release confirming PEGylation.</p>
Keywords	Polyetheylene Glycol; PEG