



# Rabbit Anti-Human SLC23A2 polyclonal antibody (CABT-L3195)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant Human Solute carrier family 23 member 2 protein (1-90AA)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	PBS with 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Long time storage is recommended at -20°C.
Ship	Wet ice

## BACKGROUND

<b>Introduction</b>	Sodium/ascorbate cotransporter. Mediates electrogenic uptake of vitamin C, with a stoichiometry of 2 Na(+) for each ascorbate.
<b>Keywords</b>	SLC23A2;solute carrier family 23 (nucleobase transporters), member 2;SLC23A1, solute carrier family 23 (nucleobase transporters), member 1;solute carrier family 23 member 2;KIAA0238;SVCT2;YSPL2;hSVCT2;yolk sac permease-like molecule 2;Na(+)/L-ascorbic acid transporter 2

## GENE INFORMATION

<b>Synonyms</b>	SLC23A2; solute carrier family 23 (nucleobase transporters), member 2; SLC23A1, solute carrier family 23 (nucleobase transporters), member 1; solute carrier family 23 member 2; KIAA0238; SVCT2; YSPL2; hSVCT2; yolk sac permease-like molecule 2; Na(+)/L-ascorbic acid transporter 2
<b>UniProt ID</b>	<a href="#">Q9UGH3</a>