



Mouse Anti-Isocarbophos monoclonal antibody, clone ICP (CABT-L4602)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The antibody recognizes the conjugated form of the small molecule. Antibody specificity was performed by using an ELISA test by competition experiments
Immunogen	Synthetic Isocarbophos conjugated to a carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ICP
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS only
Preservative	None
Storage	Store -20°C for long term.

BACKGROUND

Introduction

Isocarbophos is an organothiophosphate insecticide, an organic phosphonate, a phosphonic ester, a member of salicylates and an isopropyl ester. It has a role as an EC 3.1.1.7 (acetylcholinesterase) inhibitor, an agrochemical and an avicide. It derives from an isopropyl salicylate.

Keywords

Isocarbophos;Optunal;Isocarbofos
