



# Mouse Anti-Clodipol monoclonal antibody, clone CPD (CABT-L4605)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	The antibody recognizes the conjugated form of the small molecule. Antibody specificity was performed by using an ELISA test by competition experiments
<b>Immunogen</b>	Synthetic Clodipol conjugated to a carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	CPD
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS only
<b>Preservative</b>	None
<b>Storage</b>	Store -20°C for long term.

<b>Ship</b>	Wet ice
-------------	---------

<b>Warnings</b>	This product is for research use only and is not intended for diagnostic use.
-----------------	---

## BACKGROUND

<b>Introduction</b>	Clopidol is an organic compound that is used as in veterinary medicine as a coccidiostat. It is prepared industrially by a multistep process from dehydroacetic acid.
---------------------	---

<b>Keywords</b>	Clopidol
-----------------	----------