



Mouse Anti-Tulathromycin monoclonal antibody, clone TLM (CABT-L4609)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The antibody recognizes the conjugated form of the small molecule. Antibody specificity was performed by using an ELISA test by competition experiments
Immunogen	Synthetic Tulathromycin conjugated to a carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TLM
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS only
Preservative	None
Storage	Store -20°C for long term.

BACKGROUND

Introduction

Tulathromycin is a novel triamilide antimicrobial that has been approved for use in the treatment and prevention of bovine respiratory disease and the treatment of swine respiratory disease in the European Union and the United States.

Keywords

Tulathromycin; Tulathromycin A; C₄₁H₇₉N₃O₁₂
