



Rabbit Anti-fungal beta-glucan monoclonal antibody (CABT-RM008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	beta-glucan
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Candida (4 species), Aspergillus (4 species), P. notatum, S. cerevisiae, P. pastoris, Mucor hiemalis
Clone	9312
Purification	protein A chromatography
Conjugate	Unconjugated
Applications	IF, ELISA
Format	Liquid
Concentration	100 µg/ml
Size	100 µg
Buffer	0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide
Storage	Store at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C.

BACKGROUND

Introduction	β -Glucans (beta-glucans) comprise a group of β -D-glucose polysaccharides naturally occurring in the cell walls of cereals, bacteria, and fungi, with significantly differing physicochemical properties dependent on source.
---------------------	--

Keywords	fungus beta-glucan; fungus β -glucan; fungus
-----------------	--
