



Rabbit Anti-Streptococcus Group A carbohydrate monoclonal antibody (CABT-RM012)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	carbohydrate
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Streptococcus Group A
Clone	1546
Purification	protein A chromatography
Conjugate	Unconjugated
Applications	IF, ELISA
Format	Liquid
Concentration	100 µg/ml
Size	100 µg
Buffer	0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide
Storage	Store at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C.

BACKGROUND

Introduction

Group A streptococcus (GAS) can cause a group A streptococcal infection in hosts. Streptococcus pyogenes comprises the vast majority of the Lancefield group A streptococci, and is often used as a synonym for GAS.

Keywords

Streptococcus group A; Group A Strep; Streptococcus group A carbohydrate; Group A Strep carbohydrate
