



Rabbit Anti-Streptococcus Group B carbohydrate monoclonal antibody (CABT-RM014)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	carbohydrate
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	Streptococcus Group B
Clone	2636
Purification	protein A chromatography
Conjugate	Unconjugated
Applications	IF, ELISA
Format	Liquid
Concentration	100 µg/ml
Size	100 µg
Buffer	0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide
Storage	Store at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C.

BACKGROUND

Introduction

The bacterium *Streptococcus agalactiae* (*S. agalactiae*) is also known as group B streptococcus or GBS. Group B streptococcal infection can cause serious illness and sometimes death, especially in newborns, the elderly, and people with compromised immune systems.

Keywords

Streptococcus group B; Group B Strep; Streptococcus group B carbohydrate; Group B Strep carbohydrate
