



Mouse Anti-CDV monoclonal antibody, clone 6-5 (CABT-RM015)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	surface envelope antigen of CDV
Target	CDV
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	CDV
Clone	6-5
Purification	Chromatography on protein A Sepharose
Conjugate	Unconjugated
Applications	ELISA (Cap), PLA
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage	+4°C (2-8°C allowed)

BACKGROUND

Introduction	Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with CDV.
--------------	---

Keywords

Canine distemper virus; CDV; Canine; Canine distemper; hardpad disease
