



Mouse Anti-CPV monoclonal antibody, clone 9I8 (CABT-RM018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	CPV
Target	CPV
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	CPV
Clone	9I8
Purification	Chromatography on protein A Sepharose
Conjugate	Unconjugated
Applications	ELISA (Det), WB, ID, HIT
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage	+4°C (2-8°C allowed)

BACKGROUND

Introduction	Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunised with canine parvovirus.
--------------	---

Keywords

canine parvovirus; parvovirus; CPV; canine; parvo
