



Goat Anti-Papillomavirus Polyclonal antibody (CABT-RM239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Papillomavirus group including HPV types 11 & 16, Bovine Papillomavirus. Nonreactive with normal tissue and cultured cells
Target	Papillomavirus
Immunogen	Bovine Papillomavirus
Source/Host	Goat
Species Reactivity	BPV
Purification	>95%
Conjugate	unconjugated
Applications	ELISA, LFIA, IFA, WB, IHC
Format	Liquid
Size	500 µl
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Cold Packs

BACKGROUND

Introduction

Bovine papillomaviruses (BPV) are a paraphyletic group of DNA viruses of the subfamily Firstpapillomavirinae that are common in cattle. Infection causes warts (papillomas and fibropapillomas) of the skin and alimentary tract, and more rarely cancers of the alimentary tract and urinary bladder. They are also thought to cause the skin tumour equine sarcoid in horses and donkeys.

Keywords

Bovine papillomaviruses; BPV; papillomaviruses