



Mouse Anti-Penconazole Monoclonal Antibody, clone PCZ (CABT-Z531M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Penconazole-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PCZ
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	1-[2-(2,4-dichlorophenyl)pentyl]-1,2,4-triazole is a member of the class of triazoles that is 1,2,4-triazole substituted at position 1 by a 2-(2,4-dichlorophenyl)pentyl group. It is a dichlorobenzene and a member of triazoles.
Keywords	Penconazole; C ₁₃ H ₁₅ Cl ₂ N ₃ ; Topaze; Onmex; Topaze C; 1-[2-(2,4-dichlorophenyl)pentyl]-1,2,4-triazole