



Mouse Anti-Triazamidine Monoclonal Antibody, clone TAM (CABT-Z532M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Triazamidine-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TAM
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Triazamidine is used as a veterinary drug. Triazine has a good effect on Babesia horsebow, which can completely eliminate the worms, but it has a poor effect on Babesia horsebow.

Keywords

Triazamidine;Triazine
