



# Mouse Anti-Pyridaben Monoclonal Antibody, clone PYD (CABT-Z533M)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Pyridaben-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PYD
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

## BACKGROUND

**Introduction** Pyridaben is a pyridazinone, an organochlorine insecticide and an organochlorine acaricide. It has a role as a mitochondrial NADH:ubiquinone reductase inhibitor.

---

**Keywords** Pyridaben;96489-71-3;Sanmite

---