



Mouse Anti-Forchlorfenuron Monoclonal Antibody, clone FCF (CABT-Z534M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Forchlorfenuron-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	FCF
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice
Warnings	This product is for research use only and is not intended for diagnostic use.

BACKGROUND

Introduction	Forchlorfenuron is a plant growth regulator. It has been approved for use on kiwifruit and grapes in the United States, and it has been associated with news of watermelons exploding in China.
Keywords	Forchlorfenuron;1-(2-chloropyridin-4-yl)-3-phenylurea;1-(2-Chloro-4-pyridyl)-3-phenylurea;CPPU;N-(2-Chloro-4-pyridyl)-N'-phenylurea