



Mouse Anti-Clomazone Monoclonal Antibody, clone CMZ (CABT-Z535M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Clomazone-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CMZ
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice
Warnings	This product is for research use only and is not intended for diagnostic use.

BACKGROUND

Introduction

Clomazone is an agricultural herbicide, and has been the active ingredient of products named "Command" and "Commence". The molecule consists of a 2-chlorobenzyl group bound to a N-O heterocycle called isoxazolidinone. It is a white solid.

Keywords

Clomazone;2-[(2-chlorophenyl)methyl]-4,4-dimethyl-3-isoxazolidinone;Dimethazone
