



Mouse Anti-Oxadixyl Monoclonal Antibody, clone ODX (CABT-Z540M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Oxadixyl-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ODX
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice
Warnings	This product is for research use only and is not intended for diagnostic use.

BACKGROUND

Introduction

Oxadixyl is an oxazolidinone that is N-(2,6-dimethylphenyl)-2-methoxyacetamide in which the amide hydrogen is replaced by a 2-oxo-1,3-oxazolidin-3-yl group. A systemic fungicide used to treat seeds of a variety of food crops, as well as lawns. It has a role as an agrochemical. It is an aromatic amide, an oxazolidinone, a carbohydrazide, an ether and an anilide fungicide.

Keywords

Oxadixyl;C14H18N2O4;77732-09-3;Sandofan;Recoil;Ripost