



# Mouse Anti-Propamocarb Monoclonal Antibody, clone PPMC (CABT-Z544M)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |  |
|--------------------|--|
| Immunogen          | Propamocarb-KLH                                      |
| Isotype            | IgG  |
| Source/Host        | Mouse  |
| Species Reactivity | N/A  |
| Clone              | PPMC   |
| Purification       | Protein A purified                                   |
| Conjugate          | Unconjugated   |
| Applications       | ELISA, LFIA  |
| Format             | Liquid   |
| Concentration      | Lot specific   |
| Size               | 1 mg   |
| Buffer             | PBS  |
| Preservative       | None   |
| Storage            | Store at -20°C long term. Avoid freeze / thaw cycle. |
| Ship               | Wet ice  |

## BACKGROUND

**Introduction**

Propamocarb is a systemic fungicide used for control of soil, root and leaf disease caused by oomycetes. It is used by watering or spraying. Propamocarb is absorbed and distributed through the plant's tissue.

---

**Keywords**

Propamocarb; Propyl [3-(dimethylamino)propyl]carbamate; C<sub>9</sub>H<sub>20</sub>N<sub>2</sub>O<sub>2</sub>

---