



Mouse Anti-Flumetralin Monoclonal Antibody, clone FLM (CABT-Z545M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Flumetralin-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	FLM
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice
Warnings	This product is for research use only and is not intended for diagnostic use.

BACKGROUND

Introduction Flumetralin ia a plant growth regulator/inhibitor/modifier,it is commonly used on tobacco.

Keywords Flumetralin;C16H12ClF4N3O4
