



# Mouse Anti-Ketoconazole Monoclonal Antibody, clone KCZ (CABT-Z550M)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Ketoconazole-KLH
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	KCZ
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Ship</b>	Wet ice

## BACKGROUND

**Introduction**

Ketoconazole, sold under the brand name Nizoral among others, is an antifungal medication used to treat a number of fungal infections. Applied to the skin it is used for fungal skin infections such as tinea, cutaneous candidiasis, pityriasis versicolor, dandruff, and seborrheic dermatitis.

---

**Keywords**

Ketoconazole;Nizoral

---