



Mouse Anti-Colchicine Monoclonal Antibody, clone CCC (CABT-Z552M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Colchicine-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CCC
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Colchicine is a medication used to treat gout and Behçet's disease. In gout, it is less preferred to NSAIDs or steroids. Other uses for colchicine include the prevention of pericarditis and familial Mediterranean fever.

Keywords

Colchicine;Colcrys;Mitigare;C22H25NO6
