



Mouse Anti-Phenacetin Monoclonal Antibody, clone PNT (CABT-Z557M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Phenacetin-KLH
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PNT
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction

Phenacetin (acetophenetidin, N-(4-ethoxyphenyl)acetamide) is a pain-relieving and fever-reducing drug, which was widely used following its introduction in 1887. It was withdrawn from medicinal use as dangerous from the 1970s.

Keywords

phenacetin;Acetophenetidin;N-(4-Ethoxyphenyl)acetamide;Acetphenetidin;Acetophenetidine;Acetophenetin