



# Rabbit Anti-Human PP2D1 Polyclonal Antibody (CABT-Z662R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Protein Epitope Signature Tag (PrEST) antigen sequence: PSYQMTTDEQQIINSFYTVFREEYAAIEDLFSAINKTEAVRCEYEDTHKAFKAFWRMDR LLGLGRKEVSRVQWSGC
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Affinity Purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC Recommended dilution: IHC: 1:2500 - 1:5000
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µl
<b>Buffer</b>	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store at 4°C for short term; -20°C for long term. Avoid freeze / thaw cycle.
<b>Ship</b>	Wet ice

## Warnings

This product is for research use only and is not intended for diagnostic use.

---

# BACKGROUND

## Keywords

PP2D1;protein phosphatase 2C-like domain containing 1;C3orf48;protein phosphatase 2C-like domain-containing protein 1

---

# GENE INFORMATION

## Gene Name

PP2D1

---

## Entrez Gene ID

[151649](#)

---

## UniProt ID

[A8MPX8](#)

---