



# Anti-ID3 monoclonal antibody, clone 26F4 (CABT-ZC1175)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Monoclonal Antibody to Human ID3 molecule
<b>Antigen Description</b>	The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts.
<b>Target</b>	ID3
<b>Immunogen</b>	Full length human recombinant protein of human ID3 produced in HEK293T cell.
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	26F4
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LMNX
<b>Buffer</b>	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose.
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

## GENE INFORMATION

Entrez Gene ID [3399](#)

---

UniProt ID [Q02535](#)

---