



Anti-JUN monoclonal antibody, clone 3C3 (CABT-ZC1192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human JUN molecule
Antigen Description	JUN is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
Target	JUN
Immunogen	Full length human recombinant protein of human JUN produced in HEK293T cells.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3C3
Conjugate	Unconjugated
Applications	ELISA, LMNX
Procedure	Matched Antibody Pairs
Buffer	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose.
Preservative	0.05% Sodium Azide
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.