



Anti-NEUROG1 monoclonal antibody, clone 9H9 (CABT-ZC1239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human NEUROG1 molecule
Antigen Description	Neurogenin-1 is a protein that in humans is encoded by the NEUROG1 gene.
Target	NEUROG1
Immunogen	Full length human recombinant protein of human NEUROG1 produced in HEK293T cell.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	9H9
Conjugate	Unconjugated
Applications	ELISA, LMNX
Buffer	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol.
Preservative	0.05% Sodium Azide
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

GENE INFORMATION

Entrez Gene ID [4762](#)

UniProt ID

[E1T0H3](#)
