



# Anti-POGK monoclonal antibody, clone 4E3 (CABT-ZC1253)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Monoclonal Antibody to Human POGK molecule
<b>Antigen Description</b>	The exact function of the protein encoded by this gene is not known. However, this gene product contains a KRAB domain (which is involved in protein-protein interactions) at the N-terminus, and a transposase domain at the C-terminus, suggesting that it may belong to the family of DNA-mediated transposons in human.
<b>Target</b>	POGK
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	4E3
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, LMNX
<b>Procedure</b>	Matched Antibody Pairs
<b>Buffer</b>	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol.
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.