



Anti-WFDC2 monoclonal antibody, clone 4B0 (CABT-ZC1351)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Monoclonal Antibody to Human WFDC2 molecule
Antigen Description	This gene encodes a protein that is a member of the WFDC domain family. The WFDC domain, or WAP Signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many family members. This gene is expressed in pulmonary epithelial cells, and was also found to be expressed in some ovarian cancers. The encoded protein is a small secretory protein, which may be involved in sperm maturation.
Target	WFDC2
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	4B0
Conjugate	Unconjugated
Applications	ELISA, LMNX
Procedure	Matched Antibody Pairs
Preservative	None
Storage	Shipped at 4 °C. Upon delivery store at -20 °C. Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.