



Human BAALC blocking peptide (CDBP0549)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-BAALC antibody
Antigen Description	This gene was identified by gene expression studies in patients with acute myeloid leukemia (AML). The gene is conserved among mammals and is not found in lower organisms. Tissues that express this gene develop from the neuroectoderm. Multiple alternatively spliced transcript variants that encode different proteins have been described for this gene; however, some of the transcript variants are found only in AML cell lines.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	BAALC brain and acute leukemia, cytoplasmic [Homo sapiens]
Official Symbol	BAALC
Synonyms	BAALC; brain and acute leukemia, cytoplasmic; brain and acute leukemia cytoplasmic protein; FLJ12015;
Entrez Gene ID	79870

mRNA Refseq	NM_001024372
Protein Refseq	NP_001019543
UniProt ID	Q8WXS3
Chromosome Location	8q22.3