



Human BCAP29 blocking peptide (CDBP0570)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-BAP29 antibody
Antigen Description	BCAP29 (B-cell receptor-associated protein 29) is a protein-coding gene. Diseases associated with BCAP29 include narcolepsy, and osteoarthritis. An important paralog of this gene is BCAP31.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	BCAP29 B-cell receptor-associated protein 29 [Homo sapiens (human)]
Official Symbol	BCAP29
Synonyms	BCAP29; B-cell receptor-associated protein 29; BAP29; BCR-associated protein 29;

Entrez Gene ID	55973
mRNA Refseq	NM_001008405.2
Protein Refseq	NP_001008405.1
UniProt ID	E9PAJ1
Chromosome Location	7q22.3